

Application No. 09/908,292
Declaration Submitted Pursuant to 37 C.F.R. §1.131
Supplemental to Reply to Office Action of May 28, 2003

Agilent Dkt No. 10010218-1
R&E Dkt No. 5000-0059

BEST AVAILABLE COPY

APPENDIX (REDACTED INVENTION DISCLOSURE)



Agilent Technologies

INVENTION DISCLOSURE

PAGE ONE OF

NO 10010218 DATE RCVD

ATTORNEY MJB

Instructions: The information contained in this document is COMPANY CONFIDENTIAL and may not be disclosed to others without prior authorization. Submit this disclosure to the Agilent Technologies Legal Department as soon as possible. No patent protection is possible until a patent application is authorized, prepared, and submitted to the Government.

Descriptive Title of Invention:

Micro-fluidic flow switching device

Name of Project:

Micro Analytical Technologies

Product Name or Number: n/a

Has the invention been published, or are you planning to publish? If so, the date(s) and publication(s):

0

Has the invention been announced, offered for sale, sold, or is such activity proposed? If so, the date(s) and location(s):

0

Has the invention been disclosed to anyone outside of AGILENT TECHNOLOGIES, or will such disclosure occur? If so, the date(s) and name(s):

0

If any of the above situations will occur within 3 months, call your IP attorney or the Legal Department now at 1-553-3061 or 408-553-3061.

Has the invention been described in a lab book or other record? If so, please identify (lab book #, etc.)

088-46,47,48,49,53,54,55

Has the invention been built or tested? If so, the date:

es, [redacted]

Has this invention been made under a government contract? If so, the agency and contract number:

0

Description of Invention: Please preserve all records of the invention and attach additional pages for the following. Each additional page should be signed and dated by the inventor(s) and witness(es).

- A. Prior solutions and their disadvantages (if available, attach copies of product literature, technical articles, patents, etc.).
- B. Problems solved by the invention.
- C. Advantages of the invention over what has been done before.
- D. Description of the construction and operation of the invention (include appropriate schematic, block, & timing diagrams; drawings; samples; graphs; flowcharts; computer listings; test results; etc.)

Signature of Inventor(s): I (we) hereby submit this disclosure on this date: [redacted].

Employee No.	Name	Signature	Telex	Mailstop	Entity & Lab Name
	Kevin Kileon	[Signature]			
Employee No.	Name	Signature	Telex	Mailstop	Entity & Lab Name
	Hongfang Yin	[Signature]			
Employee No.	Name	Signature	Telex	Mailstop	Entity & Lab Name

(If more than four inventors, include additional information on another copy of this form and attach to this document)

Agilent Technologies INVENTION DISCLOSURE COMPANY CONFIDENTIAL PAGE ____ OF ____

Nature of Witness(es): (Please try to obtain the signature of the person(s) to whom invention was first disclosed.)
Invention was first explained to, and understood by, me (us) on this date: [Redacted]

Name: Tom A. van de Gort Signature: [Signature] Date of Signature: [Redacted]
Name: Harmila R. Udiavar Signature: [Signature] Date of Signature: [Redacted]

Entrepreneur & Home Address Information: (If more than four inventors, include addl. information on a copy of this form & attach to this document!)

Inventor's Full Name: [Redacted]
In Kitchen: [Redacted]
City: Palo Alto State: CA Zip: 94303
Do you have a Residential P.O. Address? P.O. BOX: No City: USA
Known as (nickname, middle name, etc.): [Redacted] Citizenship: [Redacted]

Inventor's Full Name: [Redacted]
In Kitchen: [Redacted]
City: [Redacted] State: CA Zip: 95129
Do you have a Residential P.O. Address? P.O. BOX: [Redacted]
Known as (nickname, middle name, etc.): [Redacted] Citizenship: P. R. China

Inventor's Full Name: [Redacted]
In Kitchen: [Redacted]
City: [Redacted] State: [Redacted] Zip: [Redacted]
Do you have a Residential P.O. Address? P.O. BOX: [Redacted]
Known as (nickname, middle name, etc.): [Redacted] Citizenship: [Redacted]

Inventor's Full Name: [Redacted]
In Kitchen: [Redacted]
City: [Redacted] State: [Redacted] Zip: [Redacted]
Do you have a Residential P.O. Address? P.O. BOX: [Redacted]
Known as (nickname, middle name, etc.): [Redacted] Citizenship: [Redacted]

Invention Disclosure
Micro-fluidic flow switching devices

Kevin Killeen and Hongfeng Yin, Agilent Labs, [REDACTED]

Background

Micro-fluidic (μ -fluidic) devices hold great promises for many applications. One of the most critical μ -fluidic elements is a low dead volume, reliable and leak free flow switching element (valve). Applications such as sample introduction, step gradient and linear gradient liquid chromatography (LC), the introduction of reagents into a flow stream, and the selection from different flow streams and/or fraction collection into different fluidic chambers etc. all require a valve to be successfully practiced.

Prior Art

Valve structures generated by silicon micro-machining technologies are described in US Patent Nos.: 4,768,751; 5,725,017; 5,964,239; 5,927,325; 6,102,068; 5,771,902; 5,333,831; 5,368,704; 5,417,235; 5,819,749; and 4,869,282. Laser ablated butt-coupled rotor-stator valve design has been described in US patent Nos.: 5,500,071; 5,568,413, and 5,571,410.

Summary of the invention

The present invention relates to a micro machined flow-switching device for use in micro fluidic applications. In the preferred embodiment of the present invention, two micro fluidic substrates form a high pressure capable, face seal rotary valve where fluidic structures on both devices are in fluidic communications with each other. Flow path switching is achieved with actuation and rotation of one substrate relative to the other with the capability to connect multiple flow channels sequentially or simultaneously. In another embodiment, one of the substrate is simply a rotor seal or binary shut off valve. In another embodiment, at least one of the substrates can be connected to external fluid source. An alternative embodiment would be a relative displacement (linear or combination motion e.g. slide valve), rather than just rotary motion between elements. This valve concept is extensible to multiple micro fluidic valve elements sandwiched between layers of micro fluidic structured channels.

The actuation of the valve is via mechanical coupling to a motor or pneumatically driven shaft. Alternatively, by choice of proper materials, the valve material itself may be driven by magnetic or electrical forces for non-contact actuation.

The materials of construction used in the present embodiment are low cost polymers, specifically commercial grade polyimide (Kapton™) and PEEK (polyetherether ketone). These materials can be processed by tool machining, laser ablation and plasma etching to form μ -fluidic structures in surfaces. PEEK can also be embossed or injection molded to define fluidic features. Elements constructed from alternative polymer and/or solid materials combinations are possible if they are chemically compatible with the fluids in contact with them, can withstand the internal pressures of the application without fatigue or leaks, and are mechanically stable to wear during repetitive actuation of the valve. A reliable seal between the moving elements is enhanced by using (a) compliant material(s) with the deformation of one element relative to the other forming a leak tight structure statically and during actuation.

A commercial ten-port valve design has been used to test operation in both the flow injection analysis (FIA) mode and step gradient LC mode. In step gradient mode, one loop is loaded with sample and the other loop is loaded with high percentage organic mobile phase. By switching between the two positions, sample is first loaded during the time that the LC column is flushed with mobile phase A. Then, the valve is switched. Now, sample is injected onto the LC column and non-retaining compounds are eluted by mobile phase A while Loop B is loaded with mobile phase B. Then the valve is switched again, and mobile phase B in loop B is pushed through LC column and sample is eluted through ES spray tip into Mass Spectrometer.

Functions that can be addressed by the face seal rotary valve design:


1. High pressure (>250 bar), low dead volume, leak free, simultaneous or sequential connection of multiple flow paths.
2. Fixed loop microinjection of sample for FIA mass spec analysis (commonly practiced with a 6-port design)
2. Fix loop injection for step gradient LC separation (10-port design as described above)

3. Fix loop sample injection for true gradient LC (6-port design, also with switching of splitting point)
4. Multi-position selection valve for sample selection, fraction collection and flow diversion
5. Multiple-switching with single actuation, rotor seal can address multi-channel simultaneously or actuator can address multi-rotor seals simultaneously.
6. Shut-off valve to allow or stop flow between fluid paths connected via the rotary seal.
7. Sequential switching design between differing flow paths to generate a LC gradient on chip.
8. Multilayer, stacked valve elements sandwiched between u-fluidic layered structures.



Agilent Technologies

Agilent Confidential



LC-MS Chip

The "KEY" Chip

6

Essential Elements of the Key Prototype




pump



valve



column



detector



(ESI-MS)



HP, LDV connections



Agilent Technologies

Agilent Confidential

RECEIVED
SEP 14 2004
FAX NO. 650 327 3231

Valve to Chip Interface




1 nl, laser ablated
PEEK rotor loops




Inside face of 10
port SS stator



Agilent Technologies



Kapton Key chip
face sealed between
rotor and stator.
Tested to 250 Bar.



Flow switching using 2
position actuator


Switching by valve motion
*Mechanical actuation vs.
electrical switching for CE or CEC*

Agilent Confidential


LC Column Packing

Industry-standard Packing / LC Chemistry:

Proven LC technologies and methodologies




Reflection image



50 μ m



5 μ m bead diameter




Agilent Technologies

Agilent Confidential

RECEIVED

KEY Chip Layout



Sample loop = 11.2 nl
Loop B = 26.6 nl
Column Volume before
packing = 27.3 nl

**Integrated on-chip sample
loop, mobile phase B loop,
LC column, and flow
switching means**


**Low dead and delay volume,
fast separation and washout,
low carry-over**



Agilent Technologies

Agilent Confidential

Fast Clean-up Remove DMSO interference from Reserpine



Experiment condition:

MS: Sprite single Q MS, scan
70 – 650 amu

Sample: 1 µg/ml reserpine in 0.1%
DMSO, 10 µg injected

Mobile phase A:
10% Methanol/Water, 0.5% formic
acid

Mobile Phase B:

90% Methanol/Water, 0.5% formic
acid

Flow rate: 1 µl/min

Injection: 10 nl loop

B loop: 25 nl

Step Gradient:

0 min. inject sample, load loop B

0.2 min flush column with B, new
sample loaded

0.5 min inject new sample




Agilent Technologies

Agilent Confidential



Agilent Technologies

Agenda Confidential



Step Gradient LC-MS

Experiment condition:

MS: Sprite single Q, scan 250-900 amu

Sample: Sigma LC peptide mix

Mobile phase A:

10% Methanol/Water, 0.5% formic

acid

Mobile Phase B:

90% Methanol/Water, 0.5% formic

acid

Flow rate: 1 μl/min

Injection: 10nl loop

B loop: 25ml

Step Gradient:

0 min, inject sample, load loop B

0.5 min flush column with B, new

sample loaded


12




Agilent Technologies

Agilent Confidential

Step Gradient LC-MS Clean-up and separation of peptide mixture



Flow Switching and Sample Injection for Step Gradient LC




14



Agilent Technologies

Agilent Confidential

Flow Switching and Wash for Step Gradient LC



Agilent Technologies


Agilent Confidential




Agilent Technologies

Agilent Confidential


LC-MS Concept



1. Shaft
2. Rotor
3. Caddy
4. Chip
5. Sipper
6. MS Tip
7. Backing






Rheodyne rotor seal



Rheodyne stator face assembly

RECEIVED
SEP 14 2004



2-position valve m
Hongfeng Yin

**This Page is Inserted by IFW Indexing and Scanning
Operations and is not part of the Official Record**

BEST AVAILABLE IMAGES

Defective images within this document are accurate representations of the original documents submitted by the applicant.

Defects in the images include but are not limited to the items checked:

- ☐ **BLACK BORDERS**
- ☐ **IMAGE CUT OFF AT TOP, BOTTOM OR SIDES**
- ☐ **FADED TEXT OR DRAWING**
- ☐ **BLURRED OR ILLEGIBLE TEXT OR DRAWING**
- ☐ **SKEWED/SLANTED IMAGES**
- ☐ **COLOR OR BLACK AND WHITE PHOTOGRAPHS**
- ☐ **GRAY SCALE DOCUMENTS**
- ☐ **LINES OR MARKS ON ORIGINAL DOCUMENT**
- ☐ **REFERENCE(S) OR EXHIBIT(S) SUBMITTED ARE POOR QUALITY**
- ☐ **OTHER:** _____

IMAGES ARE BEST AVAILABLE COPY.

As rescanning these documents will not correct the image problems checked, please do not report these problems to the IFW Image Problem Mailbox.